**1,25-Dihydroxyvitamin D
Cartridge Pack**

For use in conjunction with the ¹²⁵I RIA Kit

Instruction Manual

Manuel d'Instructions

Testanleitung

Manual de Instrucciones

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Bruksanvisning

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REF: 65101E

DiaSorin

Stillwater, Minnesota 55082-0285, U.S.A.

TABLE OF CONTENTS

English	Page 1
Français	2
Deutsch	3
Español	4
Italiano.....	5
Svenska	6
Norsk.....	7

1,25 DIHYDROXYVITAMIN D CARTRIDGE PACK

1. Intended Use

FOR IN VITRO DIAGNOSTIC USE.

The C₁₈OH cartridges are intended for use in conjunction with the DiaSorin 1,25-Dihydroxyvitamin D ¹²⁵I RIA kit (Cat. No. 65100E).

2. REAGENTS

Cat. No. 65101E: C₁₈OH cartridges

Upon receipt, store cartridges unopened at ambient temperature until the expiration dates listed on the labels. Once opened, the cartridges can be re-used up to 30 times or for 6 months, whichever is achieved first.

1. C₁₈OH Cartridge, 24 each.

3. PROCEDURE

3.1 Refer to the Product Insert contained within the 1,25-Dihydroxyvitamin D ¹²⁵I RIA kit (Cat. No. 65100E) for assay procedure and cartridge usage instructions.

3.2 The C₁₈OH cartridges should be washed with 5 mL of 90:10 (hexane: methylene chloride), 5 mL of IPA and 5 mL of methanol prior to first time usage. The 90:10 (hexane: methylene chloride) solvent mixture can be prepared by measuring 900 mL of hexane and 100 mL of methylene chloride. Place the new, unused cartridges on the VAC ELUT apparatus, set in the "WASTE" position, and add 5 mL of 90:10 (hexane: methylene chloride) followed by 5 mL of IPA to each column and then 5 mL of methanol. Allow each solvent mixture to pass entirely through the cartridges before proceeding to the next mixture. After this initial preparation, it will not be necessary to repeat this step since the cartridges will be regenerated during the extraction procedure.

IMPORTANT: IPA = 2-propanol or isopropanol. Do not use 1-propanol.

The C₁₈OH Cartridge material is held in place using a porous fiber frit. Discard cartridges with dislodged or missing frits as the C₁₈OH material may be lost.

4. RESULTS

1. Refer to the Product Insert contained within the 1,25-Dihydroxyvitamin D ¹²⁵I RIA kit (Cat. No. 65100E) for Results.

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.